



Laboratory Fume Hood/Canopy Plan

Florida Institute of Technology (Florida Tech)

Environmental Health & Safety (EH&S)

150 W University Blvd, Melbourne, FL 32901

Email: ehs@fit.edu

Website: <https://www.fit.edu/office-of-environmental-health-and-safety/>

REVISION HISTORY

Revision Number	Revision Date	Revised By	Description of Change
00	2019-10-09	Charles Cherrito	Initial plan creation and implementation.
01	2021-02-08	Selvin McLean	Document review; review against federal register for updated; weblinks added.

Contents

Purpose.....	4
Scope.....	4
Responsibilities.....	5
Definitions	6
Types of Hoods.....	7
Face Velocity	9
Periodic Surveys	10
General Safety Practices.....	11
Special Notes	12
References.....	13

PURPOSE

This plan is designed to assist all Florida Tech personnel in ensuring the safety of individuals who work in areas that have fume hoods or canopy-style units. Understanding the capabilities of fume hoods, their operating functions, purpose, and limitations is vital for fostering a safe work environment.

SCOPE

This plan applies to all Florida Tech locations that house fume hoods or canopy-style units.

RESPONSIBILITIES

Environmental Health & Safety (EHS)

Addresses safety needs of faculty, staff, and students through compliance assessments, chemical use surveys, and laboratory inspections. Develops and communicates protocols for laboratory and fume hood usage. Conducts risk assessment for lab design and accident investigations. Provides training assistance for fume hood use and chemical hygiene plans. The EHS Office will also conduct periodic fume hood surveys for airflow compliance.

Facilities Department

Responds to service requests regarding fume hood operations. Completes routine maintenance repairs of fume hoods and fume hood systems. Seeks the advice and approval from the EHS Office before installation of any fume hood.

Principal Investigator (PI)

Maintains user safety within the laboratory by monitoring that personnel are following procedures and that lab equipment are functioning properly; requests repair work orders from the Facilities Department to address malfunctioning systems. Maintains accurate and current information of the chemicals used in the lab and the [Chemical Hygiene Plan](#). Ensures their employees are properly trained in fume hood usage.

All Fume Hood Users

Follows all safety procedures specified by EHS, their supervisor and the [Chemical Hygiene Plan](#), this Plan, and all other pertinent SOP's/Safety Plans. Attends all required safety training sessions. Does not use fume hoods which have failed certification. Does not use fume hoods in a way they were not specifically designed. Reports fume hoods which are not operating properly, accidents, unhealthy, and unsafe conditions to the PI or supervisor

DEFINITIONS

Air Volume

When passing through a fume hood, is generally equal to the area of the sash opening multiplied by the average velocity desired. For example, if 100 feet per minute (fpm) is required and the hood has a sash opening of 7.5 square feet, then the hood's air volume is 750 (7.5 x 100) cubic feet per minute (CFM).

Air Foil

Aerodynamic sash opening directs airflow into hood and across work surface with minimum turbulence helping to ensure fume containment.

Baffle

The BAFFLE controls the pattern of the air moving into and through the fume hood

Face Velocity

Air travels in a hood at a certain speed or "velocity". The velocity of a fume hood is measured in the plane of the sash, and referred to as the face velocity, measured in feet per minute (fpm). The face velocity is related to the amount of air being pulled through the fume hood, or the volumetric rate. This number is typically measured in cubic feet per minute (CFM). The more air that is pulled through a fume hood's openings, the faster the air will travel.

Sash

The sash controls the area of the fume hood which is open. It protects the operator and controls hood face velocities. Glass options include tempered or laminated.

TYPES OF HOODS

Acid Digestion Hoods

These units are typically constructed of polypropylene to resist the corrosive effects of acids at high concentrations. If hydrofluoric acid is being used in the hood, the hood's glass sash should be constructed of polycarbonate which resists etching. Hood ductwork should be lined with polypropylene or coated with PTFE (Teflon).

Canopy Hoods

Horizontal enclosures having an open central duct suspended above a work bench or other area. Canopy hoods are most often used to exhaust areas that are too large to be enclosed within a fume hood. The major disadvantage with the canopy hood is that the contaminants are drawn directly past the user's breathing zone. Canopy hoods are only used for ventilation of heat, steam or nontoxic materials with low vapor density and not for storage or mixing of hazardous chemicals.

Explosion Proof (EP) Hoods

An EP hood is modified in a fashion to reduce spark potential outside the hood. This reduces the chances of igniting a flammable atmosphere in the room. They have specially designed EP electrical components, (switches, receptacles, and internal wiring) that are installed on site by a licensed electrician (not the manufacturer). If the room may have a potentially explosive/flammable atmosphere. ***NOTE: EP Hoods are NOT explosion proof!***

Fume Exhaust Duct Connections

Commonly called snorkels, elephant trunks or flex ducts, are designed to be somewhat mobile allowing the user to place it over a small area needing ventilation. However, for optimal efficiency, these connections must be placed within six (6) inches of an experiment, process, or equipment. These funnel-shaped exhausts aid in the removal of contaminated or irritating air from a point source to the outside.

Perchloric Acid Hoods

Concentrated or hot perchloric acid is highly oxidizing and extremely corrosive. In addition, the fumes can settle and form shock-sensitive crystals. Perchloric Acid Hoods have wash-down capabilities to prevent the buildup of explosive perchlorate salts within the exhaust systems.

Conventional Fume Hoods

Vent a constant amount of air over a 24-hour period. The volume of air passed is controlled by the sash; the higher it is opened, the higher the volume of air and the lower the face velocity. Another type of conventional hood uses a bypass to keep air volume and face velocity constant, regardless of sash height. Reduced flow fume hoods: like conventional fume hoods, also vent a constant volume of air. However, they vent air at volumes up to 50% lower than conventional fume hoods. This increases the safety of the user and reduces energy consumption.

Variable Air Volume Hoods

Use a feedback system of controls to constantly adjust the volume of air being adjusted based on sash height. The result is a constant face velocity across the front of the hood. There are several benefits to this type of hood, including decreased turbulence, increased user safety, and significantly decreased energy use. These hoods maintain the pre-selected face velocity proportional to the sash opening while monitoring hood conditions for safety.

Radioisotope Hoods

This fume hood is made with a coved stainless-steel liner and coved integral stainless-steel countertop that is reinforced to handle the weight of lead bricks. If research is planned with radioisotopes contact the [Florida Tech Radiation Safety Officer](mailto:ehs@fit.edu) at ehs@fit.edu.

MISC. NON-FUME HOODS

Although the below hoods are not considered “fume hoods”, they are common devices used research; and their usage and purpose are widely misunderstood. Therefore, they will be briefly discussed so that their functions can be distinguished.

Biosafety Cabinets (BSC's)

BSC's are enclosed, ventilated laboratory workspaces for safely working with materials contaminated with (or potentially contaminated with) pathogens requiring a defined biosafety level. They use HEPA filtration and specific airflow as a means of protection. There are many sub-types of BSC's, with each having specific characteristics as to what type of hazard can safely be manipulated in them.

Laminar Fume Hoods (LFH's)

LMF's are enclosures in which air flow is directed to prevent contamination of sterile materials. They are NOT suitable for work involving any hazards.

FACE VELOCITY

Face velocity of a very important factor in fume hoods. Airflow into a hood is achieved by an exhaust blower which “pulls” air from the laboratory room into and through the hood and exhaust system. This “pull” at the opening of the hood is measured as *face velocity*. Average face velocity is calculated by dividing the sash opening into one-foot squares. Velocity readings (feet per minute (FPM)) are taken in each grid area and averaged.

WHAT DO THE INDUSTRY STANDARDS SAY ABOUT FACE VELOCITY?

❖ **Occupational Safety and Health Administration (OSHA):**

“General air flow should not be turbulent and should be relatively uniform throughout the laboratory, with no high velocity or static areas; air flow into and within the hood should not be excessively turbulent; hood face velocity should be adequate. (Typically, 60-110 fpm.)”

❖ **ANSI/AIHA/ASSE Z9.5)**

American National Standards Institute (ANSI)
American Industrial Hygiene Association (AIHA)
The American Society of Safety Engineers (ASSE)

“Each hood shall maintain an average face velocity of 80-120 fpm with no face velocity measurement more than plus or minus 20% of average.”

❖ **Scientific Equipment and Furniture Association (SEFA):**

“Face velocities of laboratory fume hoods may be established on the basis of the toxicity or hazard of the materials used or the operations conducted within the fume hood. Note: Governmental codes rules and regulation may require specific face velocities. A fume hood face velocity of 100 fpm is considered acceptable in standard practice. In certain situations, face velocity of up to 125 fpm or as low as 75 fpm may be acceptable to meet required capture velocities of the fume hood.”

Aside from the external blower/motor unit pulling air, the face velocity is also dependent on sash height. The higher the sash height will result in less velocity (less FPM). Conversely, the lower the sash will result in an increase in velocity. This appropriate height is known as “Operating Height” and is approximately 18 inches from the base of the fume hood.

PERIODIC SURVEYS

As a part of due diligence and to ensure the safety of personnel, the EHS Office will conduct periodic surveys of all fume hoods at Florida Tech. These surveys will be performed annually, at minimum, and will consist of the following:

- Recording of average face velocity at a typical working height;
- Visual observation of any issues.

If a fume hood is not performing to standards, it will be placed out of service and the Facilities Department will be notified for addressing the issue. Upon the airflow being rectified, another survey will be performed to ensure the hood is operating as designed.

GENERAL SAFETY PRACTICES

➤ **Using Fume Hoods as Storage**

It is okay to keep chemicals in a fume hood that are being used during a procedure, but fume hoods are not designed for permanent chemical storage. Items placed in the hood can interfere with the air flow, causing turbulence which will allow contaminants to be drawn out of the hood into the room. Fume hoods are not flammables cabinets. They offer no protection from a fire occurring outside the hood in the laboratory.

➤ **Exhausting Vapors**

Fume hoods exhaust vapors into the atmosphere untreated, so it is not appropriate to use a fume hood for waste disposal. Evaporation is considered treatment of hazardous materials and is not allowed by the EPA without a permit.

➤ **Shutting Sashes Closed**

Fume hoods are designed for their sashes to be shut all the way. The airfoil sill at the base of the hood will still allow for air to be pulled into the fume hood even with the sash fully closed.

➤ **Too High of Face Velocity**

While it is important to have a face velocity between 80 fpm and 120 fpm, a higher face velocity will make the fume hood less efficient. High face velocities can create eddy (turbulent) currents that will allow for contaminants to be drawn out of the hood into the room and increase the worker's exposure.

➤ **Working Distance**

Work at least 6 inches into the hood from the face to minimize the potential for fumes to escape. As a useful reminder, place a strip of tape at this six-inch limit.

➤ **Reducing Draft**

Keep laboratory doors and windows closed and limit movement in front of the hood. Most laboratory ventilation systems are designed for labs to have all doors and windows closed. Open doors and windows can alter the air balance in the room and disrupt the airflow in the hood. Also, open doors and windows and traffic or movement in front of the hood can create turbulence resulting in vapors flowing out of the hood's interior.

➤ **Reduce Airflow Issues**

Reduce obstruction in the fume hood to improve its effectiveness by allowing adequate air flow across the working surface with minimum turbulence.

SPECIAL NOTES

- ❖ Re-circulating or ductless fume hoods are **not permitted** for the removal of chemical, biological, or radiological contaminants. Florida Tech does not service such hoods and EHS prohibits the use on campus.
- ❖ Fume hood systems **are not designed for use of biological pathogens** since they lack HEPA filtration. No biological hazards are allowed in fume hoods without EHS approval.

REFERENCES

Florida Tech Chemical Hygiene Plan

<https://www.fit.edu/office-of-environmental-health-and-safety/chemical-safety/>

29CFR 1910.1450: OSHA's

[Occupational Exposure to Hazardous Chemicals in Laboratories](#)

ANSI/ASHRAE 110-2016:

[Method of Testing Performance of Laboratory Fume Hoods](#)

ANSI/AIHA Z9.5

[The American National Standard for Laboratory Ventilation](#)

American Conference of Governmental Industrial Hygienists (ACGIH),

["Industrial Ventilation - A Manual of Recommended Practice"](#)

The American National Standard for Laboratory Ventilation (NFPA)

["Standard on Fire Protection for Laboratories Using Chemicals." Standard 45](#)